
Control of quantum sensors

#### Paola Cappellaro

Massachusetts Institute of Technology











# Limits and goals of qt. control


 Quantum information processing requires universal quantum control


- Control of quantum sensors has different goals and constraints
  - Less stringent requirements
  - New challenges:
     e.g. compromise between noise refocusing and external field sensing



# Control of quantum devices


 Continuous decoupling magnetometry






 Reconstruction of timedependent magnetic fields


 Isolated defect in diamond with electronic spin 1



• ODMR in a confocal microscope




• Electronic spin 1




- Single-spin state detection
- Fluorescence intensity (FI)
- Optical polarization
  - mK at room temperature


ODMR in a confocal microscope



Electronic spin 1



ODMR in a confocal microscope





- Single-spin state detection
  - Fluorescence intensity (FI)
- Optical polarization
- Precise control
  - Manipulation via resonant microwave



# Spin magnetometer

- Control + localization
  - → sensitivity & spatial resolution






A. Yacoby



# Spin magnetometer

- Control + localization
  - → sensitivity & spatial resolution







H.-C. Chang



Electronic spin 1




ODMR in a confocal microscope



- Single-spin state detection
  - Fluorescence intensity (FI)
- Optical polarization
- Precise control
  - Manipulation via resonant microwave



Electronic spin 1



0.5

- Single-spin state detection
  - Fluorescence intensity (FI)
- Optical polarization
- Precise control
  - Manipulation via resonant microwave

ODMR in a confocal microscope





- Up to ms, allows many operations
- Complex environment
- Nuclear spins
- Electronic Nitrogen spins



# Single-spin magnetometer

Detect magnetic field with Ramsey-type experiment





Shot-noise limited sensitivity (minimum resolvable field)


$$\delta B \sim C \frac{\hbar}{g \mu_B} \frac{1}{\sqrt{T_2^*}}$$
 [T Hz- $^{1/2}$ ]


– Noise limits the evolution time to  $au \lesssim T_2^*$ 

J. Taylor, P.C. et al., Nature Phys 2008

## Single-spin magnetometer

Detect magnetic field with Ramsey-type experiment





Shot-noise limited sensitivity (minimum resolvable field)

$$\delta B \sim C \frac{\hbar}{g \mu_B} \frac{1}{\sqrt{T_2}} \qquad \qquad \text{[T Hz-$^{1\!\!/}\!\!2]}$$

- Limited by dephasing time → Spin echo
- No sensing of DC-field, limited by pulse errors





#### CONTINUOUS DECOUPLING MAGNETOMETRY

A. Aiello, M. Hirose, P.C. Nature Comm. **4**, 1419 (2013) M. Hirose, A. Aiello, P.C. Phys. Rev. A. **86**, 062320 (2012)

 Detection of static magnetic fields Ramsey

(high sensitivity, short T<sub>2</sub>)



Detection of static magnetic fields
 Ramsey


R

high consitivity ch

Rabi\*

(high sensitivity, short  $T_2$ ) (long  $T_2$ , low sensitivity)





\*Fedder et al., Appl Phys B **102**, 497–502 (2011)



Detection of static magnetic fields Ramsey


(high sensitivity, short  $T_2$ )

Rabi\*

(long T<sub>2</sub>, low sensitivity)



**Example: Rotary Echo** 



I. Solomon (1958)



Detection of static magnetic fields


Ramsey

Rabi\*


(high sensitivity, short T<sub>2</sub>)

(long T<sub>2</sub>, low sensitivity)





**Example: Rotary Echo** 





Detection of static magnetic fields


Ramsey

(high sensitivity, short  $T_2$ )

Rabi\*

(long  $T_2$ , low sensitivity)







- Compromise:
  - Longer T<sub>2</sub> than Ramsey, higher sensitivity than Rabi
- Corrects for μw instability



## **Rotary Echo**


Intermediate (variable) T<sub>2</sub> and sensitivity








### Rotary Echo Magnetometry






# Sensitivity

Experimental results with 1 NV center




### Sensitivity

Higher sensitivity, robust against μw noise




#### Advantage: Robustness

- Rotary Echo corrects for driving fluctuations
  - Experiments (added noise)







- Adjustable coherence time
  - match external constraints, avoid overhead times




- Adjustable coherence time
  - match external constraints, avoid overhead times



- Adjustable coherence time
  - match external constraints, avoid overhead times




- Adjustable coherence time
  - match external constraints, avoid overhead times



# **AC Magnetometry**


 Pulsed and continuous Dynamical Decoupling allow measuring AC fields




- cDD allows more flexibility than PDD
  - Compromise between decoupling efficiency and bandwidth



#### **Rotary Echo**






 $\omega$ 


.20

.10

.05

#### **Rotary Echo**






 $\omega$ 


.20

.10

.05

#### **Rotary Echo**








 $\omega$ 


.20

.10

.05

#### **Rotary Echo**








 $\omega$ 

.20


.10

.05

**Rotary Echo** 



PDD (echo)



### Spectroscopy


 Frequency selectivity can be varied by pulse cycle and Rotary Echo angle

– Acquired phase is weighted by DD evolution:

$$\overline{B} = \frac{1}{t} \int_0^t b(t) F_{\text{DD}}(t) dt = \overline{B}_{\text{max}} W_{\text{DD}}(\omega)$$

### Spectroscopy

 Frequency selectivity can be varied by pulse cycle and Rotary Echo angle









#### **WAVEFORM RECONSTRUCTION**

#### Waveform reconstruction

Magnetic fields from biological activities






Magnetic activities in He-La cells


L. P. McGuinness et al., Nature Nanotechnology 6, 358–363 (2011)



CPMG sequences → frequency filters




CPMG sequences → frequency filters




Change pulse spacing

CPMG sequences → frequency filters



Change pulse spacing

CPMG sequences → frequency filters



- Limitations:
  - Inefficient reconstruction
  - Cannot reconstruct general time-dependence
  - Increasing pulse number, variable decoupling, ...



#### Waveform reconstruction

- Arbitrary waveform reconstruction with a complete (Walsh) transform
  - Complete set of orthonormal filters generated by pulsed control
  - Digital sampling reconstruction, systematic and efficient
  - Still achieves noise refocusing



#### **Walsh Functions**

Complete basis (with increasing precision @ higher order)

```
m=0 (Ramsey)

m=1 (Spin echo)

m=2 (CPMG-2)

m=3 (PDD-3)

m=4 (CPMG-4)

m=5


m=6

m=7 (PDD-7)

m=8 (CPMG-8)
```

#### **Walsh Functions**

• Complete basis (with increasing precision @ higher order)



CPMG (PDD) are a subset of Walsh



#### Walsh Reconstruction

Acquired phase depends on field and Walsh:

$$\varphi_m(T) = \gamma_e \int_0^T b(t) W_m(t/T) dt = \gamma_e \hat{b}_T(m)$$

Field obtained from Walsh transform

$$\tilde{b}_N(t) = \frac{1}{N} \sum_{m=0}^{N-1} \hat{b}_T(m) W_m(t/T)$$


Accurate reconstruction (minimal error)



#### Walsh vs. CPMG


Reconstruction of polychromatic, asymmetric

wave:  $b(t) = a_1 \sin(\omega_1 t) + a_2 \sin(\omega_2 t + \varphi)$ 



## **Preliminary results**


• Experimental reconstruction  $b(t)*\mathrm{Amp}$ 





# **Preliminary results**

• Experimental reconstruction  $b(t)*\mathrm{Amp}$ 



#### **Outlooks**

 General scheme, applicable to other measurements (e.g. temperature)

Noise spectroscopy for random fields

- Optimization
  - Optimal subset of Walsh, given prior information
  - Compressive sensing



#### Conclusions

 Quantum sensors introduce new challenges and opportunities in quantum control

- Devise control strategies to achieve flexible quantum sensors
- Quantum probe manipulation extract field information efficiently







Clarice Aiello

**Alex Cooper** 

**Gary Wolcowitz** 

Masashi Hirose

# Thanks!





Ashok Ajoy



**Honam Yum** 



Easwar Magesan





